

BLOOD BAG SHAKERS



- All-in-one
- Convenient scale design
- Easy-to-follow operation
- Alarm Indication
- Long rechargeable battery
- Logical and easy-to-use 2-button panel.
- Portable and light weight

Blood Bag Shaker MR 921 is compact instrument to provide smooth and gentle rocking. The blood collection monitor provides homogeneous mixing with anticoagulant but without the clot formation of blood cells during collection of blood from a donor. Preset standard volume is easy to change in blood collection monitor. Provision of pausing collection and change set volume during pause is also provided in blood collection mixer. Hi-Tech Instruments has provided its blood collection monitor with .This blood collection monitor also has audible alarm.

Hi-Tech Instruments
Office 7 & 8nd Floor Block C-1,
Taj Medical Complex
M.A. Jinnah Road
Karachi Pakistan.
Phone: 0092 21 2787496
Mobile: 03002343299
E-Mail: hitech.aslam@gmail.com

