

- Blood Bank Refrigerator
- Blood Bag Shaker
- Tube Stripper
- Plasma Extractor
- Platelet Incubator
- Platelet Agitator
- Platelet Rotator
- Blood Donor Chair
- Plasma Thawing Bath
- Economy Incubator
- Economy Oven

ROLLER MIXER

Various types of tubes and vials can be accommodated, with a maximum weight of 1Kg per pair of rollers. The heavy duty PVC rollers are mounted on eccentric drive spindles to provide a rolling motion, ensuring a thorough mixing action. Speed of the rollers is fixed. To cover a wide range of mixing requirements. The sprocket and gear drive, coupled with ball race roller mounting, will ensure smooth trouble free operation. The unit can run continuously for many hours unattended.

Model No. MR901

Operating Volt: 100-230V AC

Power Frequency: 50/60 Hz

Power Consumption: 5 Watt

Heavy Duty Gear Motor
Powder Coated Metal Frame
PVC Bacteria Resistant Rollers
Light Weight
Easy to Use



Hi-Tech Instruments
 Office 7 & 8 2nd Floor Block C-1,
 Taj Medical Complex
 M.A. Jinnah Road
 Karachi Pakistan.
 Phone: 0092 21 2787496
 Mobile: 03002343299
 E-Mail: hitech.aslam@gmail.com

