

- Blood Bank Refrigerator
- Blood Bag Shaker
- Tube Stripper
- Plasma Extractor
- Platelet Incubator
- Platelet Agitator
- Platelet Rotator
- Blood Donor Chair
- Plasma Thawing Bath
- Economy Incubator
- Economy Oven

## TUBE ROTATOR

Here we are introducing our Tube Rotators ideal for rotational mixing of micro centrifuge tubes, test tubes and blood tubes. These units can also be used with various temperature controlled lab equipments such as incubators, cold rooms and others with temperature range from 4° to 60°C. All these tube holders are easy to use and do not require any tools to remove or replace a tube.

Model No. TR 12  
 Operating Volt: 100-230V AC  
 Power Frequency: 50/60 Hz  
 Power Consumption: 1 Watt



Heavy Duty Gear Motor
Painted Metal Alloy Tube Disk
PVC Bacteria Resistant Exterior
Fixed Speed 10 -15 RPM
Light Weight
Easy to Use

Hi-Tech Instruments  
 Office 7 & 8 2<sup>nd</sup> Floor Block C-1,  
 Taj Medical Complex  
 M.A. Jinnah Road  
 Karachi Pakistan.  
 Phone: 0092 21 2787496  
 Mobile: 03002343299  
 E-Mail: hitech.aslam@gmail.com

