

OVEN \ INCUBATOR

Our Oven / Incubator ideal for all application with temperature range from ambient temperature 200 °C Using finely turned control technology.

Specification:-

Capacity: 35 liters

Temp. Readout: Microprocessor Controlled Digital Thermometer.


Stability: 0.5 °C.

Easy to clean interior, made of Stainless Steel with 2 Stainless Steel Grid.

Temperature Range: Ideal for incubating living culture a 37° C

However Temperature range is from Room Temperature to 200°C.

Voltage / Power Rating: 230 V, 50 Hz.



Hi-Tech Instruments
Office 7 & 8nd Floor Block C-1,
Taj Medical Complex
M.A. Jinnah Road
Karachi Pakistan.
Phone: 0092 21 2787496
Mobile: 03002343299
E-Mail: hitech.aslam@gmail.com