

PLASMA EXTRACTOR

PLASMA EXTRACTOR is designed to extract blood component from centrifuged bags. The spring loaded front panel applies pressure on the collected bag. Its sturdy to contain makes it durable and reliable, easy to use and portable. The unit is available in manual.

- Plasma Extractor Model # PE 950 is designed for automatic separation of blood components from the blood bag.
- White powder coated mild steel cabinet withstand corrosion, wear & tear on continues usage.
- Compression plate designed to exert uniform pressure on the blood bag.
- Attractive front panel bezel.
- Durable and reliable - Sturdy Construction
- Portable - Light weight
- Easy to use - Simple design



Model No.: PE 950



Model No.: PE 931

Specification	MODEL: PE950	MODEL: PE931
Power	220 V AC	_____
Sensor	Infra red	_____
Clamping	Solenoid Activated	_____
Pressure Plate	Acrylic	Acrylic
Alarm	Audio	_____
Control	Automatic	Manual

